




# What is Electroless Deposition?

In the **Electroless Deposition** process for the preparation of **bimetallic catalysts**, a reducible metal salt (cation or anion) is deposited onto the catalytically active sites of a monometallic primary metal surface through a controlled chemical reaction with a liquid-phase reducing agent. The reaction can be catalyzed by the pre-existing metal or undergo autocatalysis by the metal that is being deposited. The primary metal is typically supported on a substrate.



## Reducing agent trends for different metals

Reducing agent trends from  
Ohno, "Electrochemistry of electroless plating," Mat. Science Engr., **A146** (1991) 33-49.  
(based on electrochemical reduction of RA on different metallic electrodes)

Reducing Agent	Order of Catalytic Activity
HCHO at pH 12.5, 25C	Cu > Co ~ Au > Ag > Pt > Pd > Ni
BH <sub>4</sub> <sup>-</sup> at pH 12.5, 25C	Ni ~ Co ~ Pd > Pt > Ag ~ Au > Cu
DMAB at pH 7, 25C	Ni > Co > Pd > Pt ~ Au > Ag > Cu
N <sub>2</sub> H <sub>4</sub> at pH 12, 25C	Co > Ni > Pt ~ Pd > Cu > Ag > Au
H <sub>2</sub> PO <sub>2</sub> <sup>-</sup> at pH 9, 70C	Au > Pd ~ Ni > Co > Pt > Cu > Ag

The process is a combination of anodic and cathodic electrochemical partial reactions. The aqueous reducing agent is catalytically activated at the surface of the primary metal to produce an active hydrogen species that reduces the aqueous secondary metal salt on the surface of the primary metal.


The choice of the metal salt and the reducing agent provide a means for controlling the **Electroless Deposition** reaction. The reducing agent must be thermodynamically unstable yet kinetically stable in the absence of a catalyst, yet reactive enough to deposit a metal onto the pre-existing monometallic catalyst surface.

Controlling the pH of the aqueous bath is critical for synthesizing desirable **bimetallic catalysts** using **Electroless Deposition**. Most reducing agents are effective under basic pH conditions that are above the point of zero charge (PZC) for most supports. Therefore, strong electrostatic adsorption (SEA) of positively charged metal salts on the negatively charged support surface can be an issue. The careful selection of the positively or negatively charged form of the reducible metal salt prevents SEA.

## Effect of ligands on reduction potentials

Half-Reaction	E° (V)	Half-Reaction	E° (V)
<b>Gold</b>		<b>Copper</b>	
Au <sup>3+</sup> + 3e <sup>-</sup> → Au <sup>0</sup>	1.52	Cu <sup>2+</sup> + 2e <sup>-</sup> → Cu	0.34
AuCl <sub>4</sub> <sup>-</sup> + 3e <sup>-</sup> → Au <sup>0</sup> + 4Cl <sup>-</sup>	1.00	Cu <sup>+</sup> + e <sup>-</sup> → Cu	0.52
		Cu(CN) <sub>2</sub> <sup>-</sup> + e <sup>-</sup> → Cu + 2CN <sup>-</sup>	-0.44
Au <sup>+</sup> + e <sup>-</sup> → Au <sup>0</sup>	1.83	<b>Silver</b>	
Au(CN) <sub>2</sub> <sup>-</sup> + e <sup>-</sup> → Au <sup>0</sup> + 2CN <sup>-</sup>	-0.60	Ag <sup>+</sup> + e <sup>-</sup> → Ag	0.799
		Ag(CN) <sub>2</sub> <sup>-</sup> + e <sup>-</sup> → Ag + 2CN <sup>-</sup>	-0.31

Reaction	Log <sub>10</sub> of formation constant
AuCl <sub>3</sub> + Cl <sup>-</sup> → AuCl <sub>4</sub> <sup>-</sup>	6.0
AuCl <sub>2</sub> <sup>+</sup> + Cl <sup>-</sup> → AuCl <sub>3</sub>	4.6
Au <sup>+</sup> + 2CN <sup>-</sup> → Au(CN) <sub>2</sub> <sup>-</sup>	38.3
Cu <sup>+</sup> + 2CN <sup>-</sup> → Cu(CN) <sub>2</sub> <sup>-</sup>	24
Ag <sup>+</sup> + 2CN <sup>-</sup> → Ag(CN) <sub>2</sub> <sup>-</sup>	21.1



With appropriate **Electroless Deposition** baths, it is possible to prepare a series of **bimetallic catalysts** with increasing amounts of the deposited secondary metal and thus a range of activities and selectivities. By providing such control, **Electroless Deposition** provides a means for the rational design and synthesis of novel **bimetallic catalysts** with potentially novel catalytic properties.

Visit [www.strem.com](http://www.strem.com) for new product information and searchable catalog.


**Strem Chemicals, Inc.**  
7 Mulliken Way  
Newburyport, MA 01950-4098  
U.S.A.  
Tel.: (978) 499-1600  
Fax: (978) 465-3104  
Email: [info@strem.com](mailto:info@strem.com)

**Strem Chemicals, Inc.**  
15, rue de l'Atome  
Zone Industrielle  
67800 BISCHHEIM France  
Tel.: (33) 03 88 62 52 60  
Fax: (33) 03 88 62 26 81  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals, Inc.**  
Postfach 1215  
77672 KEHL  
Germany  
Telefon: 0 78 51/ 7 58 79  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals UK, Ltd.**  
An Independent Distributor of Strem Chemicals products  
Newton Hall, Town Street  
Newton, Cambridge  
England CB22 7ZE  
Tel: 0845 643 7263  
Fax: 0845 643 7362  
Email: [enquiries@strem.co.uk](mailto:enquiries@strem.co.uk)

## Performance of Bimetallic Catalysts prepared via Electroless Deposition



**Electroless Deposition** baths have been formulated for the deposition of both partial and multiple metal layers.

Silica-supported Ag–Pt, Au–Pd, and Cu–Pd compositions have been synthesized for use as catalysts for chemical transformations.


Carbon-supported Pt–Pd, Pt–Rh, Pt–Ru, and Pt–Co compositions have also been prepared for use in fuel cells.

**Electroless Deposition** has also been used to synthesize H<sub>2</sub>-selective Pd, Pd/Ag, and Pd/Cu alloy membranes for the production, separation and purification of ultrahigh purity H<sub>2</sub> for membrane reactor applications.

**Bimetallic catalysts** such as Ag–Pt/SiO<sub>2</sub> and Cu–Pd/SiO<sub>2</sub> have shown improved performance for the selective hydrogenation of 3,4-epoxy-1-butene (EpB), a highly functional olefin.



Au–Pd/SiO<sub>2</sub> catalysts with incremental coverages of Au on Pd have exhibited interesting catalytic performance in the hydrogenation of propylene, which is typically a structure-insensitive reaction. It was found, however, that elevated fractional coverages of Au on Pd resulted in enhanced turnover frequencies (TOFs).

In many cases, these **bimetallic catalysts** were prepared using an **Electroless Deposition** bath containing the reactive yet selective cyano salt of the secondary metal and N<sub>2</sub>H<sub>4</sub> as a strong reducing agent in combination with careful control of the pH and the reaction temperature.



## Electroless Deposition for Fuel Cell Catalysts

The commercial potential of **proton exchange membrane** (PEM) fuel cells has been limited in part by the high cost of the electrocatalysts necessary for the **oxygen reduction reaction** (ORR). The typical catalyst is platinum supported on electrically conductive porous carbon because it is corrosion resistant under the acidic conditions present in the fuel cell and is also exhibits a high activity toward oxygen reduction. High loadings (> 20 wt% Pt) are often required, however, which contributes to the high cost.



To reduce the Pt loadings, and thus the electrocatalyst cost, either the dispersion of Pt on the surface can be increased, or the specific activity of the deposited platinum for the ORR can be increased.

Using **Electroless Deposition**, it is possible to prepare platinum core/shell **bimetallic electrocatalysts** with increased specific ORR activity and stability in the corrosive PEM fuel cell environment. In one formulation, a base metal (Co) is at the core and is thus protected by the Pt shell. The core/shell structure enables a very high Pt dispersion, while possibly providing greater resistance to sintering, and allows for interaction between the Pt and the base metal atoms for enhanced activity.

Visit [www.strem.com](http://www.strem.com) for new product information and searchable catalog.

**Strem Chemicals, Inc.**  
 7 Mulliken Way  
 Newburyport, MA 01950-4098  
 U.S.A.  
 Tel.: (978) 499-1600  
 Fax: (978) 465-3104  
 Email: [info@strem.com](mailto:info@strem.com)

**Strem Chemicals, Inc.**  
 15, rue de l'Atome  
 Zone Industrielle  
 67800 BISCHHEIM France  
 Tel.: (33) 03 88 62 52 60  
 Fax: (33) 03 88 62 26 81  
 Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals, Inc.**  
 Postfach 1215  
 77672 KEHL  
 Germany  
 Telefon: 0 78 51/ 7 58 79  
 Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals UK, Ltd.**  
 An Independent Distributor of Strem Chemicals products  
 Newton Hall, Town Street  
 Newton, Cambridge  
 England CB22 7ZE  
 Tel: 0845 643 7263  
 Fax: 0845 643 7362  
 Email: [enquiries@strem.co.uk](mailto:enquiries@strem.co.uk)

The result: preparation of Pt–Co/C catalysts with reduced Pt loadings and higher oxygen reduction activity.

## Strem Products for Electroless Deposition

Strem offers a wide range of supported metals, reducible metal salts and reducing agents that are ideal for use in the synthesis of **bimetallic catalysts** via **Electroless Deposition**.

**Supported Metals** (see [www.strem.com](http://www.strem.com) for extensive list)

Support Metal	Carbon	Alumina	Barium oxide	Titania	Silica	Calcium oxide	Zinc oxide
Palladium	X	X	X	X	X	X	
Platinum	X	X			X	X	
Rhodium	X	X					
Ruthenium	X	X					
Cobalt					X		
Gold		X		X			X
Silver		X					

### Reducible Metal Salts

Gold		Palladium		Platinum		Silver		Perrhenate Salts	
NH <sub>4</sub> AuCl <sub>4</sub> ·XH <sub>2</sub> O	02-1000	Na <sub>2</sub> PdCl <sub>4</sub> ·3H <sub>2</sub> O	93-4613	Na <sub>2</sub> PtCl <sub>4</sub> ·XH <sub>2</sub> O	93-7810	KAg(CN) <sub>2</sub>	93-4739	NH <sub>4</sub> <sup>+</sup>	93-0249
NaAuBr <sub>4</sub> ·XH <sub>2</sub> O	79-3505	K <sub>2</sub> PdBr <sub>4</sub>	93-4611	K <sub>2</sub> PtBr <sub>4</sub>	78-1967				02-0900
HAuCl <sub>4</sub> ·XH <sub>2</sub> O	79-0500	Na <sub>2</sub> PdCl <sub>6</sub>	93-4612	Na <sub>2</sub> PtCl <sub>6</sub> ·6H <sub>2</sub> O	78-1995			K <sup>+</sup>	75-0500
KAu(CN) <sub>2</sub>	79-3000	K <sub>2</sub> PdCl <sub>6</sub>	93-4610	K <sub>2</sub> PtCl <sub>6</sub>	78-1960			Ag <sup>+</sup>	47-1700
KAuCl <sub>4</sub> ·XH <sub>2</sub> O	93-7906	(NH <sub>4</sub> ) <sub>2</sub> PdCl <sub>4</sub>	93-4602	K <sub>2</sub> PtCl <sub>4</sub>	78-1970	<b>Copper</b>		Na <sup>+</sup>	93-7508
		K <sub>2</sub> PdCl <sub>4</sub>	46-2126	(NH <sub>4</sub> ) <sub>2</sub> PtCl <sub>4</sub>	93-7802	KCu(CN) <sub>2</sub>	19-3025		
				K <sub>2</sub> Pt(NO <sub>2</sub> ) <sub>4</sub>	78-1985				

### Reducing Agents

Compound	Chemical Formula	Catalog No.
Tetrabutylammonium borohydride	(C <sub>4</sub> H <sub>9</sub> ) <sub>4</sub> NBH <sub>4</sub>	02-5000
Tetraethylammonium borohydride	(C <sub>2</sub> H <sub>5</sub> ) <sub>4</sub> NBH <sub>4</sub>	93-0573
Tetramethylammonium borohydride	(CH <sub>3</sub> ) <sub>4</sub> NBH <sub>4</sub>	93-0536
Sodium borohydride	NaBH <sub>4</sub>	93-1118
Sodium triacetoxyborohydride	Na(OOCCCH <sub>3</sub> ) <sub>3</sub> BH	11-0595
Sodium cyanoborohydride	NaBH <sub>3</sub> CN	11-2800
Sodium dihydrobis(2-methoxyethoxy)aluminate	Na(CH <sub>3</sub> OCH <sub>2</sub> CH <sub>2</sub> O) <sub>2</sub> AlH <sub>2</sub>	93-1077
Potassium borohydride	KBH <sub>4</sub>	93-1907
Potassium triethylborohydride	KB(CH <sub>2</sub> CH <sub>3</sub> ) <sub>3</sub> H	19-1970

If you do not see a supported metal, reducible metal salt or reducing agent of interest to you for use in the synthesis of **bimetallic catalysts** via **Electroless Deposition**, please contact Strem.

### References:

- J. Rebelli, A.A. Rodriguez, S. Ma, C.T. Williams, J.R. Monnier, Catal. Today, 160 (2011) 170.
- J. Rebelli, M. Detwiler, S. Ma, C.T. Williams, J.R. Monnier, J. Catal. 270 (2010) 224.
- M. Ohashi, K.D. Beard, S. Ma, D.A. Blom, J. St-Pierre, J.W. Van.Zee, and J.R. Monnier, Electrochimica Acta 55 (2010) 7376.
- K.D. Beard, D. Borrelli, A.M. Cramer, D. Blom, J.W. Van Zee, J.R. Monnier, ACS Nano 3 (2009) 2841.
- K.D. Beard, J.W. Van Zee, J.R. Monnier, Appl. Catal., B 88 (2009) 185.
- M.T. Schaal, A.Y. Metcalf, J.H. Montoya, J.P. Wilkinson, C.C. Stork, C.T. Williams, J.R. Monnier, Catal. Today 123 (2007) 142.
- M.T. Schaal, A.C. Pickerell, C.T. Williams, J.R. Monnier, J. Catal. 254 (2008) 131.

Visit [www.strem.com](http://www.strem.com) for new product information and searchable catalog.

**Strem Chemicals, Inc.**  
7 Mulliken Way  
Newburyport, MA 01950-4098  
U.S.A.  
Tel.: (978) 499-1600  
Fax: (978) 465-3104  
Email: [info@strem.com](mailto:info@strem.com)

**Strem Chemicals, Inc.**  
15, rue de l'Atome  
Zone Industrielle  
67800 BISCHHEIM France  
Tel.: (33) 03 88 62 52 60  
Fax: (33) 03 88 62 26 81  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals, Inc.**  
Postfach 1215  
77672 KEHL  
Germany  
Telefon: 0 78 51/ 7 58 79  
Email: [info.europe@strem.com](mailto:info.europe@strem.com)

**Strem Chemicals UK, Ltd.**  
An Independent Distributor of Strem Chemicals products  
Newton Hall, Town Street  
Newton, Cambridge  
England CB22 7ZE  
Tel: 0845 643 7263  
Fax: 0845 643 7362  
Email: [enquiries@strem.co.uk](mailto:enquiries@strem.co.uk)