


Metathesis Catalysts



STREM

Metathesis Catalysts



Strem Chemicals has been providing fine chemicals for research and commercial production for over fifty years. In this booklet you will find our selection of molybdenum and ruthenium metathesis catalysts as well as related kits for screening purposes. These catalysts can be used in a variety of metathesis applications including olefin metathesis, ring-closing metathesis, ring-opening polymerization and cross-metathesis.

At Strem, we also offer a wide variety of ligands, nanomaterials and CVD/ALD precursors. Most of our products are of high purity, typically at 99%, while some are as high as 99.9999% metals purity. We continually seek to provide new technologies from around the globe and add to our product line. We have licensing agreements with industry and academia, which allow easier access to these patent-protected products for our customers. We look forward to continued growth in order to best serve our customers' needs with the quality and service they can trust from Strem.

As part of our ongoing commitment to quality, we have achieved ISO 9001 certification for the Quality Management System (QMS) at our corporate headquarters in Newburyport, Massachusetts.

In addition, custom synthesis services are provided on a contract basis. For pharmaceutical applications, manufacturing is conducted under current Good Manufacturing Practices (cGMP) in FDA inspected kilo-lab suites. Complete documentation is available, including validation and stability studies. Active Drug Master Files (DMF's) are maintained in North America and Europe.

Our other booklets, which focus on applications and product classes, are available in print per request and also on our website. Below is a list of current booklet titles that are available. Please also check our Product Resources section online to find additional literature offerings, such as the Strem Chemiker, our technical publication, and product literature sheets.

- Biocatalysts
- Buchwald Ligands and Precatalysts
- Carbon-Based Nanomaterials & Elemental Forms
- Catalysts & Ligands
 - Sold in Collaboration with Takasago
 - Chiral Phosphoric Acids
 - Gold Elements & Compounds
 - Heterogeneous Catalysts
 - High Purity Chiral Reagents - Sold in Collaboration with Daicel
 - Kits
- Materials for Energy Applications
- Metal Catalysts for Organic Synthesis
- Metathesis Catalysts
- MOCVD, CVD & ALD Precursors
- MOFs and Ligands for MOF Synthesis
- Nanomaterials
- New Products
- Other Ligands
- Phosphorous Ligands and Compounds
- Photocatalysts
- PURATREM: High Purity Inorganics

Ephraim S. Honig, Ph.D., M.B.A
Chief Executive Officer



Metathesis 02/20
© 2020 Strem Chemicals, Inc.

Glossary of Terms

[α]D	Specific rotation
AAS	Atomic Absorption Standard
ACS	Conforms to American Chemical Society specifications
air sensitive	Product may chemically react with atmospheric oxygen or carbon dioxide at ambient conditions. Handle and store under an inert atmosphere of nitrogen or argon.
amp	Ampouled
b.p.	Boiling point in °C at 760mm, unless otherwise noted
d.	Density
dec.	Decomposes
elec. gr.	Electronic Grade, suitable for electronic applications
f.p.	Flash point in °F
gran.	Granular
heat sensitive	Product may chemically degrade if stored for prolonged periods of time at ambient temperatures or higher. Store at 5°C or lower.
hydrate	Unspecified water content which may vary slightly from lot to lot
hygroscopic	Product may absorb water if exposed to the atmosphere for prolonged periods of time (dependent on humidity and temperature). Handle and store under an inert atmosphere of nitrogen or argon.
light sensitive	Product may chemically degrade if exposed to light
liq.	Liquid
m.p.	Melting point in °C
moisture sensitive	Product may chemically react with water. Handle and store under an inert atmosphere of nitrogen or argon.
NMR grade	Suitable as a Nuclear Magnetic Resonance reference standard
optical grade	For optical applications
pwdr.	Powder
primary standard	Used to prepare reference standards and standardize volumetric solutions
PURATREM	Product has a minimum purity of 99.99% (metals basis)
purified	A grade higher than technical, often used where there are no official standards
P. Vol.	Pore volume
pyrophoric	Product may spontaneously ignite if exposed to air at ambient conditions
reagent	High purity material, generally used in the laboratory for detecting, measuring, examining or analyzing other substances
REO	Rare Earth Oxides. Purity of a specific rare-earth metal expressed as a percentage of total rare-earths oxides.
SA	Surface area
store cold	Product should be stored at -18°C or 4°C, unless otherwise noted (see product details)
subl.	Sublimes
superconductor grade	A high purity, analyzed grade, suitable for preparing superconductors
tech. gr.	Technical grade for general industrial use
TLC	Suitable for Thin Layer Chromatography
v.p.	Vapor pressure mm of Hg
xtl.	Crystalline

About Purity

Chemical purity	is reported after the chemical name, e.g. Ruthenium carbonyl, 99%
Metals purity	is reported in parentheses with the respective element, e.g. Gallium (III) bromide, anhydrous, granular (99.999%-Ga) PURATREM where 100% minus the metal purity is equal to the maximum allowable percentage of trace metal impurity

Metathesis Catalysts


MOLYBDENUM (Compounds)

42-0575

[2,6-bis(1-methylethyl)benzenaminato(2-)]
bis(1,1,1,3,3,3-hexafluoro-2-methyl-2-propano-
lato- κO)(2-methyl-2-phenylpropylidene)-, (T-4)
molybdenum, paraffin pellet (139220-25-0)
 $C_{30}H_{35}F_{12}MoNO_2$; FW: 763.53; paraffin pellet
air sensitive, moisture sensitive
Note: Developed by XiMo. Sold under license from
XiMo for research purposes only.




2pcs
10pcs

NEW



Technical Notes:

- General metathesis of many ordinary olefins, especially terminal olefins, and will ROMP many norbornene or substituted norbornadiene monomers to give all cis, and often isotactic, polymers.
- Useful for the "ring-closing" of dienes or the coupling of terminal olefins.
- Highly active and most commonly used in cross metathesis of aliphatic alkenes with 2-vinyl aromatics.



$X = S, O$

References:

- J. Am. Chem. Soc., 1994, 116, 3414.
- J. Am. Chem. Soc., 1993, 115, 4413.
- J. Am. Chem. Soc., 1992, 114, 5426.
- J. Am. Chem. Soc., 1992, 114, 7324.
- Tetrahedron, 1998, 54, 4413. (review article)
- J. Mol. Catal. Chem., 2002, 190, 45.


42-0530

[2,6-bis(1-methylethyl)benzenaminato(2-)]
[(1R)-3,3'-dibromo-2'-(1,1-dimethylethyl)
dimethylsilyl]oxy]-5,5',6,6',7,7',8,8'-octahy-
dro[1,1'-binaphthalen]-2-olate- κO](2,5-dimethyl-
1H-pyrrol-1-yl)(2-methyl-2-phenylpropylidene)
molybdenum (VI) (1103220-99-0)
 $C_{54}H_{70}Br_2MoN_2O_2Si$; FW: 1063; orange pwdr.
air sensitive, moisture sensitive
Note: Developed by XiMo. Sold under license from
XiMo for research purposes only.
Patents: U.S. 9,687,834, EP2242578.

100mg
500mg


NEW

amp




Technical Notes:

- Catalyst used in Z-selective cross metathesis.¹
- Catalyst used in enantio- or diastereo-selective ene-yne ring closing metathesis.^{1,3}



MOLYBDENUM (Compounds)

42-0530 [2,6-bis(1-methylethyl)benzenaminato(2-)][(1R)-3,3'-dibromo-2'-([(1,1-dimethylethyl)dimethylsilyl]oxy]-5,5',6,6',7,7',8,8'-octahydro[1,1'-binaphthalen]-2-olato-kO](2,5-dimethyl-1H-pyrrol-1-yl)(2-methyl-2-phenylpropylidene) molybdenum (VI) (1103220-99-0)



References:

1. *Nature*, **2008**, *456*, 933-937.
2. *J. Am. Chem. Soc.*, **2011**, *133*, 11512
3. *J. Am. Chem. Soc.*, **2009**, *131*, 10652.
4. *Org. Process Res. & Dev.*, **2016**, *20*, 1709.

42-0535 [2,6-bis(1-methylethyl)benzenaminato(2-)][(1R)-3,3'-dibromo-2'-([(1,1-dimethylethyl)dimethylsilyl]oxy)-5,5',6,6',7,7',8,8'-octahydro[1,1'-binaphthalen]-2-olato-kO](2,5-dimethyl-1H-pyrrol-1-yl)(2-methyl-2-phenylpropylidene) molybdenum (VI), paraffin pellets (75980-60-8)

NEW

$C_{54}H_{70}Br_2MoN_2O_2Si$; FW: 1063; paraffin pellet

air sensitive, moisture sensitive, (store cold)

Note: Developed by XiMo. Sold under license from XiMo for research purposes only. Patents: U.S. 9,687,834, EP2242578.

2pcs
10pcs

Technical Note:

1. See 42-0530 (page 1)

42-0520 [2,6-bis(1-methylethyl)benzenaminato(2-)]
(2,5-dimethyl-1H-pyrrol-1-yl)(2-methyl-2-phenylpropylidene)[(1,1':3',1"-terphenyl]-2'-olato)molybdenum (VI) (1703808-70-1)

NEW


amp

$C_{46}H_{50}MoN_2O$; FW: 742.86; yellow pwdr.

air sensitive, moisture sensitive, (store cold)


Note: Developed by XiMo. Sold under license from XiMo for research purposes only.

Patents: US20140309466, WO14139679.



Technical Note:

1. Catalyst for general metathesis reactions including cross metathesis and endo-selective ene-yne ring closing metathesis.¹



References:

1. *J. Am. Chem. Soc.*, **2009**, *131*, 10652.

MOLYBDENUM (Compounds)

42-0525	[2,6-bis(1-methylethyl)benzenaminato(2-)](2,5-dimethyl-1H-pyrrol-1-yl)(2-methyl-2-phenylpropylidene)[(1,1':3',1"-terphenyl]-2'-olato] molybdenum(VI), paraffin formulated pellet (1703808-70-1) C ₄₆ H ₅₀ MoN ₂ O; FW: 742.86; paraffin pellet air sensitive, moisture sensitive, (store cold) Note: Developed by XiMo. Sold under license from XiMo for research purposes only. Patents: US20140309466, WO14139679.	2pcs 10pcs
---------	---	---------------



Technical Note:

1. See 42-0520 (page 2)

42-0510	[2,6-Bis(1-methylethyl)benzenaminato(2-)(2,5-dimethyl-1H-pyrrol-1-yl)(4',6'-diphenyl[1,1':3',1"-terphenyl]-2'-olato)(2-methyl-2-phenylpropylidene)molybdenum(VI) (1572180-69-8) C ₅₆ H ₅₈ MoN ₂ O; FW: 895.05; orange powdr. air sensitive, moisture sensitive, (store cold) Note: Developed by XiMo. Sold under license from XiMo for research purposes only. Patents: US20140309466, WO14139679.	100mg 500mg
---------	---	----------------

Technical Note:

1. Catalyst for general metathesis
- ²
- reactions including cross metathesis
- ²
- and ring closing metathesis.
- ¹

Tech. Note (1)
Ref. (1)

References:

1. J. Am. Chem. Soc., 2009, 131, 10652.
2. Org. Process Res. & Dev., 2016, 20, 1709.

42-0515	[2,6-Bis(1-methylethyl)benzenaminato(2-)][(2,5-dimethyl-1H-pyrrol-1-yl)(4',6'-diphenyl[1,1':3',1"-terphenyl]-2'-olato)(2-methyl-2-phenylpropylidene)molybdenum(VI) in paraffin formulated pellet C ₅₆ H ₅₈ MoN ₂ O; FW: 895.05; paraffin pellet air sensitive, moisture sensitive, (store cold) Note: Developed by XiMo. Sold under license from XiMo for research purposes only. Patents: US20140309466, WO14139679.	2pcs 10pcs
---------	--	---------------


Technical Note:

1. See 42-0510 (page 3)

42-0560	(S)-1-((3,3'-dibromo-2'-((tert-butyldimethylsilyl)oxy)-5,5',6,6',7,7',8,8'-octahydro-[1,1'-binaphthalen]-2-yl)oxy)-N-(2,6-diisopropylphenyl)-1-(2,5-dimethyl-1H-pyrrol-1-yl)-1-(2-methyl-2-phenylpropylidene)molybdenum (VI) (1196674-83-5) C ₅₄ H ₇₀ Br ₂ MoN ₂ O ₂ Si; FW: 1063; red pwdr. air sensitive, moisture sensitive, (store cold) Note: Developed by XiMo. Sold under license from XiMo for research purposes only. Patents: U.S. 9,687,834, EP2242578.	100mg 500mg
---------	---	----------------

Technical Notes:

1. Z-selective cross metathesis.²
2. Enantio- or diastereoselective en-yne ring closing metathesis.^{1,3}



References:

1. Nature, 2008, 456, 933-937.
2. J. Am. Chem. Soc., 2011, 133, 11512.
3. J. Am. Chem. Soc., 2009, 131, 10652.

Metathesis Catalysts


MOLYBDENUM (Compounds)

42-0565 NEW	(S)-1-((3,3'-dibromo-2'-(tert-butyldimethylsilyl)oxy)-5,5',6,6',7,7',8,8'-octahydro-[1,1'-binaphthalen]-2-yl)oxy)-N-(2,6-diisopropylphenyl)-1-(2,5-dimethyl-1H-pyrrol-1-yl)-1-(2-methyl-2-phenylpropylidene)molybdenum (VI) paraffin pellets (1196674-83-5) C ₅₀ H ₆₂ Br ₂ MoN ₂ O ₂ Si; FW: 1063; paraffin pellet air sensitive, moisture sensitive, (store cold) Note: Developed by XiMo. Sold under license from XiMo for research purposes only. Patents: U.S. 9,687,834, EP2242578.	2pcs 10pcs
-----------------------	---	---------------

Technical Note:


- See 42-0560 (page 3)

42-0540 NEW amp	(R)-1-((3,3'-dibromo-2'-(tert-butyldimethylsilyl)oxy)-5,5',6,6',7,7',8,8'-octahydro-[1,1'-binaphthalen]-2-yl)oxy)-1-(2,5-dimethyl-1H-pyrrol-1-yl)-N-(2,6-dimethylphenyl)-1-(2-methyl-2-phenylpropylidene)molybdenum (VI) (1300026-28-1) C ₅₀ H ₆₂ Br ₂ MoN ₂ O ₂ Si; FW: 1006.89; orange powdr. air sensitive, moisture sensitive, (store cold) Note: Developed by XiMo. Sold under license from XiMo for research purposes only. Patents: U.S. 9,687,834, EP2242578.	100mg 500mg
------------------------------	--	----------------



Technical Note:

- Catalyst used for:
 - Enantio-selective ring opening²/cross metathesis.
 - Ring closing metathesis.
 - Enantio- and endo-selective ene-yne ring closing metathesis.



References:


- Nature, 2011, 471, 461-466.
- J. Am. Chem. Soc., 2012, 134, 2788-2799.

42-0545 NEW	(R)-1-((3,3'-dibromo-2'-(tert-butyldimethylsilyl)oxy)-5,5',6,6',7,7',8,8'-octahydro-[1,1'-binaphthalen]-2-yl)oxy)-1-(2,5-dimethyl-1H-pyrrol-1-yl)-N-(2,6-dimethylphenyl)-1-(2-methyl-2-phenylpropylidene)molybdenum (VI) paraffin pellets (1300026-28-1) C ₅₀ H ₆₂ Br ₂ MoN ₂ O ₂ Si; FW: 1006.89; paraffin pellet air sensitive, moisture sensitive, (store cold) Note: Developed by XiMo. Sold under license from XiMo for research purposes only. Patents: U.S. 9,687,834, EP2242578.	2pcs 10pcs
-----------------------	--	---------------

Technical Note:

- See 42-0540 (page 4)

42-1205 amp	2,6-Diisopropylphenylimidoneophylidene molybdenum(VI) bis(hexafluoro-t-butoxide) SCHROCK'S CATALYST (139220-25-0) Mo(C ₁₀ H ₁₂)(C ₁₂ H ₁₇ N)(OC(CH ₃)(CF ₃) ₂) ₂ ; FW: 765.53; yellow to orange powdr. air sensitive, moisture sensitive, (store cold)	100mg 500mg 2g
----------------	--	----------------------



Technical Notes:




- Unlike Mo(C₁₀H₁₂)(C₁₂H₁₇N)(OC₄H₉)₂, the bis(hexafluoro-t-butoxide) (MoF₆) derivative will metathesize many ordinary olefins, especially terminal olefins, and will ROMP many norbornene or substituted norbornadiene monomers to give all cis, and often isotactic, polymers.

Metathesis Catalysts

MOLYBDENUM (Compounds)

42-1205 2,6-Diisopropylphenylimidoneophylidene molybdenum(VI) bis(hexafluoro-t-butoxide)
(continued) SCHROCK'S CATALYST (139220-25-0)

2. Useful for the "ring-closing" of dienes or the coupling of terminal olefins.
3. Useful for cross-metathesis of aliphatic alkenes with 2-vinyl aromatics.




X = CH = CH, S, O

References:



1. J. Am. Chem. Soc., **1994**, 116, 3414.
2. J. Am. Chem. Soc., **1993**, 115, 4413.
3. J. Am. Chem. Soc., **1992**, 114, 5426.
4. J. Am. Chem. Soc., **1992**, 114, 7324.
5. Tetrahedron, **1998**, 54, 4413. (review article)
6. J. Mol. Catal. Chem., **2002**, 190, 45.

42-0550 (2,5-Dimethyl-1H-pyrrol-1-yl)
NEW
amp (2,2",4,4",6,6"-hexamethyl[1,1':3',1"-terphenyl]-
2'-olato)(2-methyl-2-phenylpropylidene)
[2,3,4,5,6-pentafluorobenzenaminat(2-κN]
,(T-4) molybdenum(VI) (1433803-79-2)
C₄₉H₄₅F₆MoN₂O; FW: 832.81; orange pwdr.
air sensitive, moisture sensitive, (store cold)
Note: Developed by XiMo. Sold under license from
XiMo for research purposes only.
Patent: U.S. 9,441,059.



Technical Notes:

1. Z-selective alkene halide synthesis by cross metathesis.¹
2. Stereoselective macrocyclic ring closing metathesis.²



Macrocyclic Z-enolates

References:

1. Nature, **2016**, 531, 459.
2. J. Am. Chem. Soc., **2014**, 136, 16493.

MOLYBDENUM (Compounds)

42-0555 NEW	(2,5-Dimethyl-1H-pyrrol-1-yl)(2,2'',4,4'',6,6''-hexamethyl[1,1':3',1''-terphenyl]-2'-olato)(2-methyl-2-phenylpropylidene)[2,3,4,5,6-pentafluorobenzenaminato(2-)-kN](T-4) molybdenum(VI), paraffin pellets (1433803-79-2) $C_{46}H_{44}F_5MoN_2O$; FW: 832.81; paraffin pellet air sensitive, moisture sensitive, (store cold) Note: Developed by XiMo. Sold under license from XiMo for research purposes only. Patent: U.S. 9,441,059.	2pcs 10pcs
----------------	--	---------------

Technical Note:



- See 42-0550 (page 5)

RUTHENIUM (Compounds)

96-0440	Apeiron Ammonium Catalysts Kit See page 23	
96-0420	Apeiron Bulky Catalysts Kit See page 24	
96-0410	Apeiron nitro-Grela Catalysts Kit See page 25	
96-0430	Apeiron Polymerization Catalysts Kit See page 26	
96-0400	Apeiron Ruthenium Metathesis Catalyst Kit See page 27	
44-0775 NEW	Bis(1-(2,6-diethylphenyl)-3,5,5-trimethyl-3-phenylpyrrolidin-2-ylidene)(3-phenyl-1H-inden-1-ylidene)ruthenium(II) dichloride UltraCat (2055540-61-7) $C_{61}H_{68}Cl_2N_2Ru$; FW: 1001.18; burgundy pwdr. air sensitive, (store cold) Note: Sold in collaboration with Apeiron Synthesis, Inc. Patent: PCT/IB2016/054486.	100mg 500mg

Technical Note:

- Highly selective versatile cross metathesis catalyst leading to terminal and internal olefins, such as ethenolysis via ring closing metathesis (including challenging macrocyclization), cross metathesis (involving electron deficient partners like acrylates), or ene-yne metathesis. These transformations are highly selective and take place at exceptionally low catalyst loadings.




References:

- Angew. Chem. Int. Ed. Engl., 2017, 56, 981-986.

44-0759	(1,3-Bis(2,6-diisopropylphenyl)-4-((4-ethyl-4-methylpiperazin-1-ium-1-yl)methyl)imidazolidin-2-ylidene)(2-isopropoxybenzylidene)ruthenium(II) dichloride chloride dihydrate FixCat (1799947-97-9) $C_{45}H_{66}Cl_2N_4ORuCl\cdot 2(H_2O)$; FW: 887.47 (923.50); green pwdr. (store cold) Note: Sold in collaboration with Apeiron Synthesis, Inc. U.S. Patent 61/603,790; PCT/EP2013/053967	100mg 500mg
---------	--	----------------

Technical Note:


- Fixcat is a stable olefin metathesis initiator with very good solubility in neat water. The product efficiently promotes ring-closing, cross, and enyne metathesis reactions of water soluble substrates. Suitable for homogeneous and heterogeneous in batch or flow setup.



RUTHENIUM (Compounds)

44-0759 (1,3-Bis(2,6-diisopropylphenyl)-4-((4-ethyl-4-methylpiperzain-1-ium-1-yl)methyl)imidazolidin-2-ylidene)(2-isopropoxybenzylidene)ruthenium(II)dichloride chloride dihydrate FixCat (1799947-97-9)


The Fixcat structure was also optimized to serve as a versatile and very stable catalyst, easily immobilized on solid supports. In its SCA-15 supported version, it showed exceptional efficiency in promoting ring-closing and cross-metathesis reactions, in both batch and continuous flow setups. Fix Cat is also applicable as a homogeneous catalyst, where compatible solvents include alcohols and halogenated solvents.



References:


1. *ChemSus Chem.*, **2015**, *8*, 4139.

44-0055 [1,3-Bis(2,6-di-i-propylphenyl)-4,5-dihydroimidazol-2-ylidene]-[2-i-propoxy-5-(trifluoroacetamido)phenyl]methyleneruthenium(II) dichloride M71-S1Pr (1212008-99-5)
 $C_{39}H_{50}Cl_2F_3N_3O_2Ru$; FW: 821.80; green pwdr.
Note: Sold under license from Omega Cat System for research purposes only
WO 2008/065187, PCT/EP2008/054901, Fr n°08/05403.



Technical Note:


1. Enhanced activity ruthenium "boomerang" pre-catalyst used in the olefin metathesis, enyne metathesis, and cross metathesis reactions, that can be recycled, and leaves reduced ruthenium in the product after silica gel chromatography.



References:


1. *Eur. J. Org. Chem.*, **2009**, 4254.
2. *J. Org. Chem.*, **2008**, *73*, 4225.

44-0750 [1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene][2-[(1-(methoxy(methyl)amino)-1-oxopropan-2-yl)oxy]benzylidene]ruthenium(II) dichloride GreenCat (1448663-06-6)
 $C_{39}H_{52}Cl_2N_2O_3Ru$; FW: 783.33; green pwdr.
Note: Sold in collaboration with Aperion Synthesis, Inc. U.S. Patent 61/666,009, PCT/EP2013/062435. Aperion Ruthenium Metathesis Catalyst Kit component.




Technical Note:

1. Efficient, durable and reusable olefin metathesis catalyst with a high affinity to silica gel.



RUTHENIUM (Compounds)

44-0750 [1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene]{2-[1-(methoxy(methyl)amino)-1-oxopropan-2-yl]oxy}benzylidene)ruthenium(II) dichloride GreenCat (1448663-06-6)
 (continued)




Tech. Note (1)
 Ref. (1)

References:

1. *Tetrahedron*, 2013, 69, 7408

44-0748 [1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene][(2-[(1-methoxy(methyl)amino)-1-oxopropan-2-yl]oxy]benzylidenediodoruthenium(II) GreenCat-I2
 C₃₉H₅₃I₂N₃O₃Ru; FW: 966.74; green solid
air sensitive, (store cold)
 Note: Recommended storage 2-8°C. Sold in collaboration with Apeiron Synthesis, Inc.
 U.S. Patent 6/1666,009 PCT/EP2013/062435.




Technical Note:

1. Catalyst for metathesis applications.

References:

1. *Angew. Chem.* 2017, 129, 1001 –1006.

44-0770 1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene)(2-i-propoxy-5-nitrobenzylidene) ruthenium(II) dichloride Nitro-Grela SiPr (928795-51-1)
 C₃₇H₄₈Cl₂N₃O₃Ru; FW: 755.78; green pwdr.
air sensitive, (store cold)
 Note: Recommended storage 2-8°C. Sold in collaboration with Apeiron Synthesis, Inc.
 U.S. Patent 6/867,303 PCT/EP2003/0112




Technical Note:

1. Catalyst for metathesis applications.

References:

1. *ACS Catalysis*, 2017, 7(8), 5443-5449.
2. *Catalysis Science & Technology* 2017, 7(6), 1284-1296.
3. *Angew. Chem.* 2017, 129, 1001 –1006.

44-0782 [1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene)(2-i-propoxy-5-nitrobenzylidene) ruthenium(II) diiodide nitro-Grela I2 SiPr (1874265-00-5)
 C₃₇H₄₉I₂N₃O₃Ru; FW: 938.68; olive brown pwdr.
air sensitive, (store cold)
 Note: Recommended storage 2-8°C. Sold in collaboration with Apeiron Synthesis, Inc.
 U.S. Patent 6/867,303 PCT/EP2003/01122.



Technical Note:

1. Catalyst for metathesis applications.


References:

1. *Beilstein Journal of Organic Chemistry*, 2015, 11, 1823-1832.

RUTHENIUM (Compounds)

44-0793

[1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene]tricyclohexylphosphine)-(2-oxobenzylidene)ruthenium(II) chloride LatMet S1Pr (1544328-59-7)
 $C_{59}H_{72}ClN_2OPRu$; FW: 913.68; dark green xtl.
air sensitive, (store cold)
 Note: Recommended storage 2-8°C. Sold in collaboration with Apeiron Synthesis, Inc.
 U.S. Patent 9,328,132, PCT/EP2013/065839.

100mg
500mg

NEW

Technical Note:


1. Catalyst for metathesis applications.

References:

1. *U.S. Pat. Appl. Publ. 2015*, US 20150158896 A1 20150611.
2. *PCT Int. Appl. 2014*, WO 2014016422 A1 20140130.




44-0944

Benzylidene-bis(tricyclohexylphosphine)dichlororuthenium, 97%
 [Grubbs Catalyst Gen. 1] (172222-30-9)
 $C_{43}H_{72}Cl_2P_2Ru$; FW: 822.97; purple pwdr.; m.p. 153 (dec)
(store cold)

1g
5g
25g

Technical Notes:

1. N-Heterocyclic carbene and phosphine ruthenium indenylidene precatalysts: A comparative study in olefin metathesis.
2. Ruthenium(IV) alkylidenes as precatalysts for direct arylations of alkenes with aryl chlorides and an application to sequential catalysis.
3. Synthesis of high molecular weight comb block copolymers and their assembly into ordered morphologies in the solid state.
4. Autonomic healing of epoxy vinyl esters via ring opening metathesis polymerization.
5. Amide synthesis from alcohols and amines catalyzed by ruthenium N-heterocyclic carbene complexes.
6. Stereoselective synthesis of E,Z-configured 1,3-dienes by ring-closing metathesis. Application to the total synthesis of lactimidomycin.
7. A highly efficient olefin metathesis process for the synthesis of terminal alkenes from fatty acid esters.
8. Long-Spaced aliphatic polyesters.
9. Enaminones via rhodium-catalyzed coupling of thioamides and α -diazocarbonyl compounds.
10. Stereoselective synthesis of macrocyclic peptides via a dual olefin metathesis and ethenolysis approach.
11. Synthesis and characterization of precisely-defined ethylene-co-aryl ether polymers via ADMET polymerization.
12. Regio- and stereoselective dehydrogenative silylation and hydrosilylation of vinylarenes catalyzed by ruthenium alkylidenes.


Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)Tech. Note (3)
Ref. (3)

Metathesis Catalysts


RUTHENIUM (Compounds)

**44-0944
(continued)**


Benzylidene-bis(tricyclohexylphosphine)dichlororuthenium, 97% [Grubbs Catalyst Gen. 1]
(172222-30-9)




Tech. Note (4)
Ref. (4)




Tech. Note (5)
Ref. (5)




Tech. Note (7)
Ref. (7)




Tech. Note (8)
Ref. (8)




Tech. Note (9)
Ref. (9)



Tech. Note (10)
Ref. (10)



Tech. Note (11)
Ref. (11)



Tech. Note (12)
Ref. (12)


RUTHENIUM (Compounds)

44-0944 Benzylidene-bis(tricyclohexylphosphine)dichlororuthenium, 97% [Grubbs Catalyst Gen. 1]
(continued) (172222-30-9)

References:



1. *Chem. Eur. J.*, **2007**, *13*, 8029.
2. *Angew. Chem. Int. Ed.*, **2007**, *46*, 6364.
3. *J. Am. Chem. Soc.*, **2007**, *129*, 10551.
4. *Adv. Funct. Mater.*, **2008**, *18*, 44.
5. *Chem. Eur. J.*, **2010**, *16*, 6820.
6. *J. Am. Chem. Soc.*, **2011**, *133*, 9232.
7. *Top. Catal.*, **2012**, *55*, 518.
8. *Macromolecules*, **2013**, *46*, 7213.
9. *J. Org. Chem.*, **2014**, *79*, 7405.
10. *Chem. Sci.*, **2015**, *6*, 4561
11. *Polymer*, **2015**, *64*, 76.
12. *Org. Lett.*, **2016**, *18*, 5324.

44-0073 Bis(tricyclohexylphosphine)[(phenylthio)methylene]ruthenium(II) dichloride, min. 97% (219770-99-7)
 $\text{RuCl}_2(\text{CHSC}_6\text{H}_5)_2[\text{P}(\text{C}_6\text{H}_{11})_3]_2$; FW: 855.02;
purple pwdr.

1g
5g

Technical Notes:


1. Metathesis catalyst, stable in air and can be used in aqueous media.
2. Catalyst of choice for the ring-opening metathesis polymerization of cycloolefins.
3. Catalyst concentration 2-3 times lower than comparable phenyl and vinyl substituted ruthenium carbenes.
4. Excellent initiator for solvent-free polymerization and control of initiation rates and gelation times.
5. Highly selective catalyst for the ring opening/cross-metathesis of norbornene derivatives.

Tech. Note (5)
Ref. (3)

References:


1. *J. Organomet. Chem.*, **2000**, *606*, 16
2. "Ring Opening Metathesis Polymerization of Related Chemistry", NATO Science Series II, **2002**, 23.
3. *Organometallics*, **2003**, *22*, 586

44-0063 Bis(tricyclohexylphosphine)-3-phenyl-1*H*-inden-1-ylideneruthenium(II) dichloride (250220-36-1)
 $\text{RuCl}_2(\text{C}_15\text{H}_{10})[\text{P}(\text{C}_6\text{H}_{11})_3]_2$; FW: 923.07; brown pwdr.
Note: Sold in collaboration with Umicore for research purposes only.

1g
5g

Technical Notes:

1. Highly active, air-stable catalyst used for the ring closing metathesis of dienes.
2. Used in cross-metathesis.

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (2)

References:


1. *Angew. Chem. Int. Ed.*, **2000**, *112*, 3012. (review article)
2. *Adv. Synth. Catal.*, **2008**, *350*, 2959.

RUTHENIUM (Compounds)

44-0755

NEW

1,3-Bis(2,4,6-trimethylphenyl)-4-[(trimethylammonio)methyl]imidazolidin-2-ylidene]-2-(i-propoxybenzylidene)dichlororuthenium(II) hexafluorophosphate StickyCat PF6
 $C_{35}H_{46}Cl_2F_3N_3OPRu$; FW: 843.72; green pwdr.
air sensitive, (store cold)
 Note: Recommended storage 2-8°C. Sold in collaboration with Apeiron Synthesis, Inc.
 U.S. Patent 61/603,790, PCT/EP2013/053967.

100mg
500mg

Technical Note:

1. Catalyst for metathesis applications.


Reference:

1. *Angew. Chem.* **2017**, 129, 1001 –1006.

44-0795

NEW

1,3-Bis(2,4,6-trimethylphenyl)-4-[(trimethylammonio)methyl]imidazolidin-2-ylidene]-2-(i-propoxy-5-nitrobenzylidene)dichlororuthenium(II) chloride nitro-StickyCat Cl (1415661-45-8)
 $C_{35}H_{47}Cl_3N_4O_3Ru$; FW: 779.20; green pwdr.
air sensitive, (store cold)
 Note: Recommended storage 2-8°C. Sold in collaboration with Apeiron Synthesis, Inc.
 U.S. Patent 61/603,790 PCT/EP2013/053967.

100mg
500mg

Technical Note:


1. Catalyst for olefin metathesis applications.

References:

1. *U.S. Pat. Appl. Publ.* **2013**, US 20130225807 A1 20130829.
2. *Green Chemistry*, **2012**, 14(12), 3264-3268.


44-0083

1,3-Bis(2,4,6-trimethylphenyl)-4,5-dihydroimidazol-2-ylidene[2-(i-propoxy)-5-(N,N-dimethylaminosulfonyl)phenyl]methyleneruthenium(II) dichloride (resin supported) Zhan Catalyst II
 FW: >1000; black solid; Loading: 0.5 mmol/g
 Note: Sold under license from Zannan for research purposes only. Patents CN1907992A, US 2007/0043180 A1, PCT WO 2007/003135 A1.

100mg
500mg
2g

Technical Notes:

1. Efficient, air-stable metathesis catalyst.
2. Used in cross-metathesis/Friedel-Crafts Cascade reaction.


Tech. Note (2)
Ref. (1)

References:

1. *Angew. Chem. Int. Ed.*, **2009**, 48, 7428.

44-0082

1,3-Bis(2,4,6-trimethylphenyl)-4,5-dihydroimidazol-2-ylidene[2-(i-propoxy)-5-(N,N-dimethylaminosulfonyl)phenyl]methyleneruthenium (II) dichloride, Zhan Catalyst-1B, min 96% (918870-76-5)
 $RuCl_2[C_{21}H_{26}N_2][C_{12}H_{17}NO_3S]$; FW: 733.75; green solid
 Note: Sold under license from Zannan for research purposes only. Patents CN1907992A, US 2007/0043180 A1, PCT WO 2007/003135 A1.

100mg
500mg
2g


Technical Notes:

1. See 44-0083 (page 12)

RUTHENIUM (Compounds)

44-0768 [1,3-Bis(2,4,6-trimethylphenyl)-4-[(4-ethyl-4-methylpiperazin-1-ium-1-yl)methyl]imidazolidin-2-ylidene]-{2-i-propoxybiphenylidene}dichlororuthenium(II) chloride AquaMet (1414707-08-6)
 $C_{39}H_{55}Cl_3N_3ORu$; FW: 803.31; green pwdr.


Note: Sold in collaboration with Apeiron Synthesis, Inc.
 U.S. Patent 61/603,790, PCT/EP2013/053967.
 Aperion Ruthenium Metathesis Catalyst Kit component.




100mg
 500mg

Technical Notes:

1. This is an highly active catalyst for olefin metathesis in water.
2. This catalyst is used as linker free heterogeneous metathesis catalysts after immobilization on silica or siliceous mesoporous molecular sieves. TON of ~16000 for the RCM of (-)- β -citronellene are reported.²



Tech. Note (1)
 Ref. (1)




Tech. Note (1)
 Ref. (1)

References:

1. *Catal. Sci. Technol.*, **2012**, 22, 2424
2. *ACS Catal.*, **2014**, 4, 3227


44-0047 [1,3-Bis(2,4,6-trimethylphenyl)-2-imidazolidinylidene]-[2-[[4-methylphenyl]imino]methyl]-4-nitrophenolyl]-[3-phenyl-1H-inden-1-ylidene]ruthenium(II) chloride (934538-04-2)
 $C_{50}H_{47}ClN_4O_3Ru$; FW: 888.46; orange-brown solid
 Note: Sold in collaboration with Umicore for research purposes only.
 EP 1 468 004 B1, US 2002/0349956.




100mg
 500mg

Technical Note:

1. Highly active catalyst for ring opening polymerization and ring closing metathesis reactions.



Ref. (1-3)



Ref. (1-3)

References:

1. WO 2003062253.
2. *Eur. J. Org. Chem.*, **2009**, 655.
3. *Curr. Org. Synth.*, **2008**, 5, 291


RUTHENIUM (Compounds)

44-0049

[1,3-Bis(2,4,6-trimethylphenyl)-2-imidazolidinylidene]-[2-[(2-methylphenyl)imino] methyl]phenolyl]-[3-phenyl-1H-inden-1-ylidene]ruthenium(II) chloride
(934538-12-2)


$C_{50}H_{48}ClN_3ORu$; FW: 843.46; red-brown solid
Note: Sold in collaboration with Umicore for research purposes only.

EP 1 468 004 B1, US 2002/0349956.


100mg
500mg

Technical Note:

- Highly active catalyst for ring opening and ring closing metathesis reactions.



Ref. (1-3)



Ref. (1-3)

References:

- WO 2003062253.
- Eur. J. Org. Chem., 2009, 655.
- Curr. Org. Synth., 2008, 5, 291


44-0026

1,3-Bis(2,4,6-trimethylphenyl)-2-imidazolidinylidene)(3-phenyl-1H-inden-1-ylidene)(4,5-dichloro-1,3-diethyl-1,3-dihydro-2H-imidazol-2-ylidene)ruthenium(II) dichloride (1228169-92-3)

$C_{43}H_{46}Cl_4N_4Ru$; FW: 861.73; orange-brown solid



Note: Sold in collaboration with Umicore for research purpose only.

Patent US 10,873,026.

50mg
250mg
1g


Technical Note:

- Catalyst used for RCM reactions leading to trisubstituted olefins.

Tech. Note (1)
Ref. (1)Tech. Note (1)
Ref. (1)

RUTHENIUM (Compounds)

44-0026 1,3-Bis(2,4,6-trimethylphenyl)-2-imidazolidinylidene)(3-phenyl-1H-inden-1-ylidene)(4,5-dichloro-1,3-diethyl-1,3-dihydro-2H-imidazol-2-ylidene)ruthenium(II) dichloride
 (continued) (1228169-92-3)




Tech. Note (1)
 Ref. (1)

References:

1. *Organometallics*, **2010**, *29*, 2761

44-0767 [1,3-Bis(2,4,6-trimethylphenyl)imidazolidin-2-ylidene]-2-i-propoxy-5-nitrobenzylidene)ruthenium(II) diiodide nitro-Grela I2
 (1874264-99-9)
 $C_{31}H_{37}I_2N_3O_3Ru$; FW: 854.52; olive brown pwdr.
air sensitive, (store cold)
 Note: Recommended storage 2-8°C. Sold in collaboration with Apeiron Synthesis, Inc.
 U.S. Patent 6/867,303 PCT/EP2003/01122.




Technical Note:

1. Catalyst for metathesis applications.

References:


1. *Beilstein Journal of Organic Chemistry* **2015**, *11*, 1823-1832.

44-0758 [1,3-Bis(2,4,6-trimethylphenylimidazolidin-2-ylidene)]-2-i-propoxy-5-nitrobenzylidene)ruthenium(II) dichloride nitro-Grela (502964-52-5)
 $C_{31}H_{37}Cl_2N_3O_3Ru$; FW: 671.62; green pwdr.
 Note: Sold in collaboration with Apeiron Synthesis, Inc. U.S. Patent 6,867,303, PCT/EP2003/01122. Aperion Ruthenium Metathesis Catalyst Kit component.



Technical Notes:

1. This catalyst is used for an efficient synthesis of HCV Protease Inhibitor, BILN 2061.
2. The nitro-Grela catalyst is a stable olefin metathesis initiator active in various ring-closing, cross, and enyne metathesis reactions. Its efficiency has been proven in numerous total syntheses of natural and biologically active compounds, and in material science (see Ref. 1-9)



Tech. Note (1)
 Ref. (1,2)


References:

1. *Org. Lett.*, **2008**, *10*, 1303
2. *Org. Process Res. Dev.*, **2009**, *13*, 250
3. *Org. Lett.*, **2013**, *15*, 1016
4. *Chem. Eur. J.*, **2013**, *19*, 11847
5. *Org. Lett.*, **2010**, *12*, 248
6. *Angew. Chem. Int. Ed.*, **2008**, *47*, 6483
7. *Tetrahedron Lett.*, **2006**, *47*, 6351
8. *J. Am. Chem. Soc.*, **2004**, *126*, 9318
9. *Angew. Chem. Int. Ed.*, **2002**, *41*, 4038

RUTHENIUM (Compounds)


44-0753

[1,3-Bis(2,4,6-trimethylphenyl)imidazolidin-2-ylidene](tricyclohexylphosphine)-(2-oxobenzylidene)ruthenium(II) chloride LatMet (1407229-58-6)
 $C_{46}H_{66}ClN_2OPRu$; FW: 829.52; green pwdr.
 Note: Sold in collaboration with Apeiron Synthesis, Inc. U.S. Patent 9,328,132, PCT/EP2013/065839. Aperion Ruthenium Metathesis Catalyst Kit component.

100mg
500mg

Technical Note:

1. This catalyst is used for olefin metathesis polymerization.

Tech. Note (1)
Ref. (1,2)


References:

1. *Chem. Eur. J.*, **2014**, *20*, 14120
2. *Eur. J. Inorg. Chem.*, **2014**, 1131

44-0787

NEW

[1,3-Bis(2,4,6-trimethylphenyl)imidazolidin-2-ylidene](tricyclohexylphosphine)-(2-oxo-5-nitrobenzylidene)ruthenium(II) chloride Nitro-LatMet (1544328-53-1)
 $C_{46}H_{64}ClN_3O_3OPRu$; FW: 874.52; brown xtl.
air sensitive, (store cold)
 Note: Recommended storage 2-8°C. Sold in collaboration with Apeiron Synthesis, Inc. U.S. Patent 9,328,132, PCT/EP2013/065839.

100mg
500mg

Technical Note:


1. Catalyst for metathesis applications.

References:

1. *Angew. Chem.* **2017**, *129*, 1001–1006.
2. *U.S. Pat. Appl. Publ.* **2015**, US 20150158896 A1 20150611.
3. *Chemistry - A European Journal*, **2014**, *20*(43), 14120-14125.
4. *PCT Int. Appl.* **2014**, WO 2014016422 A1 20140130.
5. *European Journal of Inorganic Chemistry*, **2014**, *2014*(7), 1131-1136.


44-0765

[1,3-Bis(2,4,6-trimethylphenyl)-4-[(trimethylammonio)methyl]imidazolidin-2-ylidene]-2-(i-propoxybenzylidene)dichlororuthenium(II) chloride StickyCat Cl (1452227-72-3)
 $C_{35}H_{48}Cl_3N_3ORu$; FW: 734.20; green pwdr.
 Note: Sold in collaboration with Apeiron Synthesis, Inc. U.S. Patent 61/603,790, PCT/EP2013/053967. Aperion Ruthenium Metathesis Catalyst Kit component.

100mg
500mg

Technical Note:

1. This is an easily removable olefin metathesis catalyst.
 It shows high activity at 40-110°C. Residual ruthenium metal is usually below 5 ppm.

Tech. Note (1)
Ref. (1)

References:

1. *Green Chem.*, **2012**, *14*, 3264

RUTHENIUM (Compounds)

44-0797 NEW	Dichloro[1,3-Bis(2,6-di-i-propylphenyl)-4-((4-ethyl-4-methylpiperazin-1-ium-1-yl)methyl)imidazolidin-2-ylidene](2-isopropoxybenzylidene)ruthenium(II) hexafluorophosphate FixCat PF6 $C_{45}H_{67}Cl_2F_6N_4OPRu$; FW: 996.98; green pwdr. air sensitive, (store cold) Note: Recommended storage 2-8°C. Sold in collaboration with Apeiron Synthesis, Inc. U.S. Patent 61/603,790 PCT/EP2013/053967.		100mg 500mg
-----------------------	---	--	----------------

Technical Note:

1. Catalyst for metathesis applications.


References:

1. *Angew. Chem.* **2017**, *129*, 1001–1006.

44-0760	Dichloro[1,3-bis(2,4,6-trimethylphenyl)imidazolidin-2-ylidene]2-[(ethoxy-2-oxoethylidene)amino]benzylidene)ruthenium(II) HeatMet $C_{33}H_{37}Cl_2N_3O_2Ru$; FW: 667.63; dark purple xtls. (store cold) Note: Sold in collaboration with Apeiron Synthesis, Inc. U.S. Patent 14/443,034; PCT/IN2013/002543.		100mg 500mg
---------	--	--	----------------

Technical Note:

1. HeatMet catalyst is a highly efficient latent catalyst, requiring thermal activation to initiate catalytic activity. Its characteristics are especially suitable to mold polymerization of reactive monomers such as dicyclopentadiene (DCPD). The catalyst is soluble in toluene and dichloromethane.

Tech. Note (1)
Ref. (1)

References:

1. *ChemSus Chem.*, **2015**, *8*, 4139.

44-0792 NEW	Dichloro[1,3-di-i-propylphenylimidazolidin-2-ylidene]2-[(ethoxy-2-oxoethylidene)amino]benzylidene)ruthenium(II) HeatMet SIPr (2097273-88-4) $C_{35}H_{49}Cl_2N_3O_2Ru$; FW: 751.79; dark violet pwdr. air sensitive, (store cold) Note: Recommended storage 2-8°C. Sold in collaboration with Apeiron Synthesis, Inc. U.S. Patent 14/443,034.		100mg 500mg
-----------------------	--	--	----------------

Technical Note:

1. Catalyst for metathesis applications.

References:

1. *ACS Catalysis*, **2017**, *7*(6), 4115–4121.


44-0778 NEW	(1-(2,6-Diethylphenyl)-3,5,5-trimethyl-3-phenylpyrrolidin-2-ylidene)(2-isopropoxy-5-nitrobenzylidene)ruthenium(II) dichloride UltraNitroCat (2106819-64-9) $C_{33}H_{40}Cl_2N_2O_3Ru$; FW: 684.66; green pwdr. air sensitive, (store cold) Note: Sold in collaboration with Apeiron Synthesis, Inc. Patent: PCT/IB2016/054486.		100mg 500mg
-----------------------	---	--	----------------

Technical Note:

1. High turnover selective catalyst for cross metathesis in macrocyclization reactions.

RUTHENIUM (Compounds)

44-0778 (1-(2,6-Diethylphenyl)-3,5,5-trimethyl-3-phenylpyrrolidin-2-ylidene)(2-isopropoxy-5-nitrobenzylidene)ruthenium(II) dichloride UltraNitroCat (2106819-64-9)
 (continued)



Tech. Note (1)
Ref. (1)

GC yield = 90%; TON = 30000


References:

1. ACS Catal., 2017, 7, 5443-5449

44-0740 (1,3-Di-o-tolylimidazolidin-2-ylidene)(2-i-propoxy-5-nitrobenzylidene)dichlororuthenium(II)

NEW Nitro-Grela Si-o-Toly

C₂₇H₂₅Cl₂N₃O₃Ru; FW: 615.51; green pwdr.
air sensitive, (store cold)
Note: Recommended storage 2-8°C. Sold in
collaboration with Apeiron Synthesis, Inc.
U.S. Patent 6/867,303 PCT/EP2003/01122.



100mg
500mg

Technical Note:


1. Catalyst for metathesis applications.

References:

1. Angew. Chem. 2017, 129, 1001 –1006.

44-7778 3-Phenyl-1H-inden-1-ylidene[bis(i-butylphobane)] ruthenium(II) dichloride (894423-99-5)


C₃₉H₅₆Cl₂P₂Ru; FW: 758.78; red pwdr.
Note: Sold in collaboration with Umicore for research
purposes only. Patent US 10,518,716.



250mg
1g

Technical Note:

1. This catalyst exhibits high selectivity as a general purpose metathesis catalyst for applications other than polymerization. It has improved air, moisture and heat resistance.




References:

1. J. Org. Chem., 2008, 73, 259

44-0078 {[2-(i-Propoxy)-5-(N,N-dimethylaminosulfonyl)phenyl]methylene}(tricyclohexylphosphine)ruthenium(II) dichloride Zhan Catalyst -1C (918871-44-0)


C₃₀H₅₆Cl₂NO₂PRuS; FW: 707.74; brown solid
Note: Sold under license from Zannan for research
purposes only. Patents CN1907992A,
US 2007/0043180 A1, PCT WO 2007/003135 A1.



500mg
2g

Technical Notes:

1. Efficient, air-stable metathesis catalyst.
2. Used in cross-metathesis/Friedel-Crafts Cascade reaction.




Tech. Note (1)

Metathesis Catalysts

RUTHENIUM (Compounds)


44-0078 {2-(i-Propoxy)-5-(N,N-dimethylaminosulfonyl)phenyl]methylene}(tricyclohexylphosphine)
 (continued) ruthenium(II) dichloride Zhan Catalyst -1C (918871-44-0)



References:



1. PCT Int. Appl. (2007), WO 2007003135 A1
2. *Angew. Chem. Int. Ed.*, **2009**, *48*, 7428

44-7785 Tricyclohexylphosphine[1,3-bis(2,4,6-trimethylphenyl) imidazol-2-ylidene][2-thienylmethylene]ruthenium(II) dichloride, min. 95% [catMETium® RF 2] (1190427-49-6)
 $C_{44}H_{61}Cl_2N_2PRuS$; FW: 852.98;
 brown pwdr.
 Note: Sold in collaboration with Evonik for research purposes only.
 Patent US 6635768.



Technical Note:


1. Efficient catalyst for ring-closing metathesis.



References:

1. *Organometallics*, **1999**, *18*, 5416
2. *Chemistry Today*, **2009**, *27*, 24

44-7780 Tricyclohexylphosphine[1,3-bis(2,4,6-trimethylphenyl)-4,5-dihydroimidazol-2-ylidene][(phenylthio)methylene]ruthenium(II) dichloride (1155422-69-7)
 $RuCl_2[C_{21}H_{26}N_2][C_6H_5S][P(C_6H_{11})_3]$; FW: 881.04; purple-brown solid





Technical Notes:

1. Metathesis catalyst, stable in air and can be used in aqueous media.
2. Catalyst of choice for the ring-opening metathesis polymerization of cycloolefins.
3. Catalyst concentration 2-3 times lower than comparable phenyl and vinyl substituted ruthenium carbenes.
4. Excellent initiator for solvent-free polymerization and control of initiation rates and gelation times.
5. Highly selective catalyst for the ring opening/cross-metathesis of norbornene derivatives.

Metathesis Catalysts


RUTHENIUM (Compounds)

44-7780 Tricyclohexylphosphine[1,3-bis(2,4,6-trimethylphenyl)-4,5-dihydroimidazol-2-ylidene] [(phenylthio)methylene]ruthenium(II) dichloride (1155422-69-7)





References:

1. *J. Org. Chem.*, **2000**, *60*, 65.
2. "Ring Opening Metathesis Polymerization of Related Chemistry", NATO Science Series II, **2002**, 23.
3. *Organometallics*, **2003**, *22*, 586.




Technical Note:

1. Efficient catalyst for ring-closing metathesis.



References:

1. *Organometallics*, **1999**, *18*, 5416
2. *Chemistry Today*, **2009**, *27*, 24



Technical Note:

1. See 44-7775 (page 20)

RUTHENIUM (Compounds)

44-7790	Tricyclohexylphosphine[4,5-dimethyl-1,3-bis(2,4,6-trimethylphenyl)imidazol-2-ylidene][2-thienylmethylene]ruthenium(II) dichloride, min. 95% [catMETium® RF 3] (1190427-50-9) C ₄₆ H ₆₅ Cl ₂ N ₂ PRuS; FW: 881.04; violet to brown pwdr. Note: Sold in collaboration with Evonik for research purposes only. Patent US 6635768.		100mg 500mg 2g
---------	---	--	----------------------


Technical Note:

- See 44-7775 (page 20)

44-7777	Tricyclohexylphosphine[3-phenyl-1H-inden-1-ylidene][1,3-bis(2,4,6-trimethylphenyl)-4,5-dihydroimidazol-2-ylidene]ruthenium(II) dichloride, min. 95% (536724-67-1) C ₅₄ H ₆₉ Cl ₂ N ₂ PRu; FW: 949.09; dark red pwdr. Note: Sold in collaboration with Umicore for research purposes only. For use in lifescience applications and research purposes only.		100mg 500mg
---------	---	--	----------------

Technical Note:

- Catalyst used for the Atom Transfer Radical Polymerization (ATRP) of vinyl monomers.



References:

- New J. Chem., 2003, 27, 257

44-0763	Tricyclohexylphosphine(2-i-propoxy-5-nitrobenzylidene)dichlororuthenium(II) Nitro-Grela 1 gen. (625082-83-9) C ₂₈ H ₄₄ Cl ₂ NO ₃ PRu; FW: 645.60; brown pwdr. air sensitive, (store cold) Note: Recommended storage 2-8°C. Sold in collaboration with Apeiron Synthesis, Inc. U.S. Patent 6/867,303 PCT/EP2003/01122.		100mg 500mg
---------	--	--	----------------

Technical Note:

- Catalyst for metathesis applications.


References:

- ACS Catalysis, 2017, 7(8), 5443-5449.
- Angew. Chem. 2017, 129, 1001 –1006.
- Chemistry - A European Journal, 2014, 20(42), 13716-13721.
- Advanced Synthesis & Catalysis, 2013, 355(10), 1997-2006.

RUTHENIUM (Compounds)




44-7783

Tri(i-propoxy)phosphine(3-phenyl-1H-inden-1-ylidene)[1,3-bis(2,4,6-trimethylphenyl)-4,5-dihydroimidazol-2-ylidene]ruthenium (II) dichloride, min. 95%
cis-Caz-1 (1255536-61-8)
 $C_{45}H_{51}Cl_2N_2O_3PRu$; FW: 876.89; brown pwdr.

50mg
250mg

Technical Notes:

1. Solvent dependent catalyst for the photoinduced ring opening metathesis polymerization of 1,5 cyclooctadiene
2. Thermo and Photoactivation of 44-7783 at elevated and room temperatures enables a wide variety of olefin metathesis reaction to be carried out under mild conditions with excellent yields

Tech. Note (1)
Ref. (1)Tech. Note (2)
Ref. (1,2)

References:

1. ACS Catal., 2018, 8, 6413–6418.
2. Chem. Commun., 2010, 7115.

Metathesis Catalysts

KITS - Apeiron Ammonium Catalysts Kit

96-0440

Apeiron Ammonium Catalysts Kit

NEW

Sold in collaboration with Apeiron Synthesis, Inc. Complexes from our Ammonium Catalysts Kit can be applied to metathesis in neat water [44-0768 (AquaMet), 44-0765 (StickyCat Cl), 44-0795 (nitro-StickyCatCl)], ethyl acetate/dimethyl carbonate [44-0755 (StickyCat PF6), 44-0797 (FixCat PF6)] or heterogeneously after catalyst deposition on solid support [44-0768 (AquaMet), 44-0759 (FixCat)].

These catalysts are especially recommended for applications in which low levels of residual ruthenium is desired.

Components also available for individual sale. Contains the following:

 44-0755	 44-0759	 44-0765
100mg	100mg	100mg
 44-0768	 44-0795	 44-0797
100mg	100mg	100mg

44-0755	1,3-Bis(2,4,6-trimethylphenyl)-4-[(trimethylammonio)methyl]imidazolidin-2-ylidene]-[2-i-propoxybenzylidene] dichlororuthenium(II) hexafluorophosphate StickyCat PF6	100mg	See page 12
44-0759	(1,3-Bis(2,6-diisopropylphenyl)-4-((4-ethyl-4-methylpiperazin-1-ium-1-yl)methyl)imidazolidin-2-ylidene)-(2-isopropoxybenzylidene)ruthenium(II)dichloride chloride dihydrate FixCat (1799947-97-9)	100mg	See page 6
44-0765	[1,3-Bis(2,4,6-trimethylphenyl)-4-[(trimethylammonio)methyl]imidazolidin-2-ylidene]-[2-i-propoxybenzylidene] dichlororuthenium(II) chloride StickyCat Cl (1452227-72-3)	100mg	See page 16
44-0768	[1,3-Bis(2,4,6-trimethylphenyl)-4-[(4-ethyl-4-methylpiperazin-1-ium-1-yl)methyl]imidazolidin-2-ylidene]-[2-i-propoxybenzylidene] dichlororuthenium(II) chloride AquaMet (1414707-08-6)	100mg	See page 13
44-0795	1,3-Bis(2,4,6-trimethylphenyl)-4-[(trimethylammonio)methyl]imidazolidin-2-ylidene)-(2-i-propoxy-5-nitrobenzylidene)dichlororuthenium(II) chloride nitro-StickyCat Cl (1415661-45-8)	100mg	See page 12
44-0797	Dichloro(1,3-Bis(2,6-di-i-propylphenyl)-4-((4-ethyl-4-methylpiperazin-1-ium-1-yl)methyl)imidazolidin-2-ylidene)-(2-isopropoxybenzylidene)ruthenium(II) hexafluorophosphate FixCat PF6	100mg	See page 17

KITS - Apeiron Bulky Catalysts Kit


96-0420

Apeiron Bulky Catalysts Kit

NEW

Sold in collaboration with Apeiron Synthesis, Inc. Catalysts within the Bulky Catalysts Kit were designed to be less sensitive to minor impurities that are commonly present in metathesis substrates/solvents. The bulkiness of the ligands in this kit helps to reduce the risk of unwanted double-bond migration. These complexes are especially recommended for RCM and CM of sterically non-demanding substrates.

Components also available for individual sale. Contains the following:



44-0748	[1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene][(2-(1-methoxy(methyl)amino)-1-oxopropan-2-yl)oxy]benzylidenediiodoruthenium(II) GreenCat-I2	100mg	See page 8
44-0750	[1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene]{2-[(1-methoxy(methyl)amino)-1-oxopropan-2-yl]oxy}benzylidene)ruthenium(II) dichloride GreenCat (1448663-06-6)	100mg	See page 7
44-0770	1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene)(2-i-propoxy-5-nitrobenzylidene)ruthenium(II) dichloride Nitro-Grela SiPr (928795-51-1)	100mg	See page 8
44-0782	[1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene](2-i-propoxy-5-nitrobenzylidene)ruthenium(II) diiodide nitro-Grela I2 SiPr (1874265-00-5)	100mg	See page 8

Metathesis Catalysts

KITS - Apeiron nitro-Grela Catalysts Kit













96-0410

Apeiron nitro-Grela Catalysts Kit

NEW

Sold in collaboration with Apeiron Synthesis, Inc. Our nitro-Grela catalysts kit contains complexes with wide-ranging activity and application profiles. These testing catalysts have excellent potential within the early stages of development.

Components also available for individual sale. Contains the following:

 44-0740	 44-0758	 44-0763
 44-0767	 44-0770	 44-0782
 44-0740	 44-0758	 44-0763
 44-0767	 44-0770	 44-0782

44-0740	(1,3-Di-o-tolylimidazolidin-2-ylidene)(2-i-propoxy-5-nitrobenzylidene)dichlororuthenium(II) Nitro-Grela Si-O-Tolyl	100mg	See page 18
44-0758	[1,3-Bis(2,4,6-trimethylphenylimidazolidin-2-ylidene)]-(2-i-propoxy-5-nitrobenzylidene)ruthenium(II) dichloride nitro-Grela (502964-52-5)	100mg	See page 15
44-0763	Tricyclohexylphosphine(2-i-propoxy-5-nitrobenzylidene)dichlororuthenium(II) Nitro-Grela 1 gen. (625082-83-9)	100mg	See page 21
44-0767	[1,3-Bis(2,4,6-trimethylphenyl)imidazolidin-2-ylidene]-(2-i-propoxy-5-nitrobenzylidene)ruthenium(II) diiodide nitro-Grela I2 (1874264-99-9)	100mg	See page 15
44-0770	1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene)-(2-i-propoxy-5-nitrobenzylidene)ruthenium(II) dichloride Nitro-Grela SiPr (928795-51-1)	100mg	See page 8
44-0782	[1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene]-(2-i-propoxy-5-nitrobenzylidene)ruthenium(II) diiodide nitro-Grela I2 SiPr (1874265-00-5)	100mg	See page 8

Metathesis Catalysts

KITS - Apeiron Polymerization Catalysts Kit


96-0430

Apeiron Polymerization Catalysts Kit

NEW


Sold in collaboration with Apeiron Synthesis, Inc. The catalysts included in our Polymerization Catalysts Kit were specifically designed for ROMP of strained monomers such as dicyclopentadiene or norbornene. The latency of these catalysts allows for controlled preparation of formulation and curing within various conditions.

Components also available for individual sale. Contains the following:




44-0753

100mg




44-0760

100mg




44-0787

100mg



44-0792

100mg



44-0793

100mg

44-0753	[1,3-Bis(2,4,6-trimethylphenyl)imidazolidin-2-ylidene](tricyclohexylphosphine)-(2-oxobenzylidene)ruthenium(II) chloride LatMet (1407229-58-6)	100mg	See page 16
44-0760	Dichloro[1,3-bis(2,4,6-trimethylphenyl)imidazolidin-2-ylidene]{2-[(ethoxy-2-oxoethylidene)amino]benzylidene} ruthenium(II) HeatMet	100mg	See page 17
44-0787	[1,3-Bis(2,4,6-trimethylphenyl)imidazolidin-2-ylidene](tricyclohexylphosphine)-(2-oxo-5-nitrobenzylidene)ruthenium(II) chloride Nitro-LatMet (1544328-53-1)	100mg	See page 16
44-0792	Dichloro[1,3-di-i-propylphenyl]imidazolidin-2-ylidene}{2-[(ethoxy-2-oxoethylidene)amino]benzylidene} ruthenium(II) HeatMet SIPr (2097273-88-4)	100mg	See page 17
44-0793	[1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene](tricyclohexylphosphine)-(2-oxobenzylidene)ruthenium(II) chloride LatMet SIPr (1544328-59-7)	100mg	See page 9


Metathesis Catalysts

KITS - Apeiron Ruthenium Metathesis Catalyst Kit

96-0400


Apeiron Ruthenium Metathesis Catalyst Kit

Sold in collaboration with Apeiron Synthesis, Inc. This kit also includes a metal scavenger (07-2203). Components also available for individual sale. Contains the following:




07-2203

1g




44-0750

100mg




44-0753

100mg




44-0758

100mg




44-0759

100mg




44-0760

100mg



44-0765

100mg



















44-0768

100mg

07-2203	1,4-Bis(2-isocyanopropyl)piperazine (SnatchCat Metal Scavenger) (51641-96-4)	1g	Visit strem.com
44-0750	[1,3-Bis(2,6-di-i-propylphenyl)imidazolidin-2-ylidene]-{2-[1-(methoxy(methyl)amino)-1-oxopropan-2-yl]oxy}benzylidene)ruthenium(II) dichloride GreenCat (1448663-06-6)	100mg	See page 7
44-0753	[1,3-Bis(2,4,6-trimethylphenyl)imidazolidin-2-ylidene]-{2-[1-(ethoxy-2-oxoethylidene)amino]benzylidene}(tricyclohexylphosphine)-(2-isopropoxybenzylidene)ruthenium(II) dichloride LatMet (1407229-58-6)	100mg	See page 16
44-0758	[1,3-Bis(2,4,6-trimethylphenyl)imidazolidin-2-ylidene]-{2-i-propoxy-5-nitrobenzylidene)ruthenium(II) dichloride nitro-Grela (502964-52-5)	100mg	See page 15
44-0759	(1,3-Bis(2,6-diisopropylphenyl)-4-[(4-ethyl-4-methylpiperazin-1-ium-1-yl)methyl]imidazolidin-2-ylidene)(2-isopropoxybenzylidene)ruthenium(II)dichloride chloride dihydrate FixCat (1799947-97-9)	100mg	See page 6
44-0760	Dichloro(1,3-bis(2,4,6-trimethylphenyl)imidazolidin-2-ylidene)-{2-[(ethoxy-2-oxoethylidene)amino]benzylidene)ruthenium(II) HeatMet	100mg	See page 17
44-0765	[1,3-Bis(2,4,6-trimethylphenyl)-4-[(trimethylammonio)methyl]imidazolidin-2-ylidene]-{2-(i-propoxybenzylidene)}dichlororuthenium(II) chloride StickyCat Cl (1452227-72-3)	100mg	See page 16
44-0768	[1,3-Bis(2,4,6-trimethylphenyl)-4-[(4-ethyl-4-methylpiperazin-1-ium-1-yl)methyl]imidazolidin-2-ylidene]-{2-(i-propoxybenzylidene)}dichlororuthenium(II) chloride AquaMet (1414707-08-6)	100mg	See page 13

Available Booklets



Visit strem.com for new product announcements.

The Strem Product Line

OUR LINE OF RESEARCH CHEMICALS

Biocatalysts & Organocatalysts	Custom Synthesis
Electronic Grade Chemicals	cGMP Facilities
Fullerenes	FDA Inspected
High Purity Inorganics & Alkali Metals	Drug Master Files
Ionic liquids	Complete Documentation
Ligands & Chiral Ligands	
Metal Acetates & Carbonates	
Metal Alkoxides & beta-Diketonates	
Metal Alkyls & Alkylamides	
Metal Carbonyls & Derivatives	
Metal Catalysts & Chiral Catalysts	
Metal Foils, Wires, Powders & Elements	
Metal Halides, Hydrides & Deuterides	
Metal Oxides, Nitrates & Chalcogenides	
Metal Scavengers	
Metallocenes	
Nanomaterials	
Organofluorines	
Organometallics	
Organophosphines & Arsines	
Porphines & Phthalocyanines	
Precious Metal & Rare Earth Chemicals	
Volatile Precursors for MOCVD, CVD & ALD	

Strem Chemicals, Inc.

7 Mulliken Way
Dexter Industrial park
Newburyport, MA 01950-4098
U.S.A

Office Tel: (978) 499-1600
Office Fax: (978) 465-3104

Toll-free (U.S. & Canada)
Tel: (800) 647-8736
Fax: (800) 517-8736

Email: info@strem.com
www.strem.com

Strem Chemicals, Inc.

15, rue de l'Atome
Zone Industrielle
F-67800 BISCHHEIM (France)
Tel.: +33 (0) 3 88 62 52 60
Fax: +33 (0) 3 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals, Inc.

Postfach 1215
D-77672 KEHL, Germany
Tel.: +49 (0) 7851 75879
Fax: +33 (0) 3 88 62 26 81
Email: info.europe@strem.com

Strem Chemicals UK, Ltd.

An Independent Distributor
of Strem Chemicals Products
Newton Hall, Town Street
Newton, Cambridge, CB22 7ZE, UK
Tel: +44 (0)1223 873 028
Fax: +44 (0)1223 870207
Email: enquiries@strem.co.uk

